
**Mini Vortex Stirrer LBX Instruments, V03 series, without
speed control**

Please read the User Manual carefully before use, and follow all operating and safety instructions!

user manual

english / español / français / italiano

User manual

LBX Instruments V03 series Mini Vortex Stirrer, without speed control

Introduction

Users must read this manual carefully before using this item. Follow up the instructions and procedures to be aware of all the precautions, and in order to obtain the maximum benefits and equipment durability.

Service

If help is needed, contact with your supplier/manufacturer or Labbox via: **www.labbox.com**

Please, provide the customer service representative with the following information:

- Serial number (on the back side)
- Description of the problem
- Your contact information

Warranty

The Mini Vortex Mixer has a 24 months warranty from the purchase date for materials and/or production faults if used as described in this manual.

This warranty is not applied to the equipment or to any damaged part due to a bad installation, bad connections, wrong use, accident, or abnormal use conditions.

For warranty claims, please contact your supplier/manufacturer.

1. Safety instructions

	Connect the equipment to a grounded power supply to ensure the safety of the equipment and the experiment. Connect the power when the equipment requires it.
	The use of this equipment in inflammable and explosive, toxic, or highly corrosive experiments is forbidden.
	Place the Vortex mixer on a flat, horizontal, and steady surface, leaving at least 30 cm of free space on each side.
	The Vortex mixer must be only used by qualified staff aware of the equipment's handling in this manual.
	Do not overfill the testing tubes, only testing tubes with stopper can be used with this device.

- Make sure that only trained staff operates this equipment.
- In case of abnormal conditions during the mixer's performance, please, unplug immediately the device and contact your supplier/manufacturer.
- While using the equipment, use personal protective equipment according to the substance's mixing hazard category to prevent the risks from:
 - Splashing and liquids evaporation
 - Toxic or flammable gas releases
 - Projected parts
- Do not use this equipment in explosive and flammable environments, or with hazardous materials.
- The equipment can only be disconnected from the power supply by pulling out the socket, do not pull the wire.
- Make sure that the power supply's voltage matches the necessary electric input for this device.

2. Use Instructions

This equipment is designed for mixing substances with educational purposes, in laboratories, or industry. This article is not suitable for use in residential areas or dangerous environments.

3. Inspection

3.1. Unpacking

Unpack the device carefully and check it for any apparent damages that may have arisen during transportation. If necessary, please contact your supplier.

Note:

If there is any apparent damage to the equipment, do not plug it into the power line.

3.2. Items list

The package includes the following items:

Items	Qty
Main unit	1
User Manual	1
Power adapter	1

4. Device parts

1. Base plate
2. Soft area for testing tubes

Figure 1: VORX-003-001 series

5. Operating

- Place the equipment on a flat and steady surface that can support the equipment's weight leaving 30 free cm on each side.
- Connect the device's power adapter to the electric network.
- Push gently the test tube into the Vortex soft area.
- To stop the process, just lift up the test tube.

6. Faults

- The equipment does not start.
 - Make sure the power supply is working.
 - Make sure the cable is plugged in.

If these faults are not solved, contact the supplier / manufacturer.

7. Cleaning and maintenance

- The Vortex mixer is maintenance-free.
- Do not spray cleaning agents directly onto the equipment. Always use an impregnated cloth or tissue
- Disconnect the equipment from the power supply while cleaning.
- Use only the recommended cleaning agents for each situation:

Dyes	Isopropyl alcohol
Construction materials	Water containing tensides / Isopropyl alcohol
Cosmetics	Water containing tensides / Isopropyl alcohol
Food products	Water containing tenside
Fuels	Water containing tenside

- Before using any other cleaning agents or decontamination methods, verify with the manufacturer or supplier if the chosen method will not harm the equipment. Protective gloves shall be worn while cleaning.

Note:

- Electronic devices may not be cleaned with cleanser.
- If maintenance service is needed, the equipment must be cleaned prior to shipping, to avoid pollution from hazardous substances, and shipped in its original packaging.
- If the equipment will not be used for a long time, power off and place in a dry, clean, steady location at room temperature.

8. Technical Characteristics

Items	Parameters
Voltage [V]	100-240
Frequency [Hz]	50/60
Power [W]	1.2W
Speed range [rpm]	300-2000
Max. [H ₂ O,L]	1.5
Base plate [mm]	Ø120
Base plate material	PET (tereftalato de polietileno)
Container material	Anti-inflamable
Dimensions [mm] D x W x H	100 x 78
Weight (kg)	0.55
Room temperature permitted [°C]	5-40
Relative humidity permitted	80%
Protection class (DIN EN60529)	II

9. Accessories

Please contact our company to order the accessories available for this device: testing tubes with stopper

For additional help, please visit www.labbox.com .

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