V01 Mini vortex stirrer, without speed control

Please read the User Manual carefully before use, and follow all operating and safety instructions!

user manual

english

User Manual

V01 Mini vortex stirrer, without speed control

Preface

Users should read this Manual carefully, follow the instructions and procedures, and beware of all the cautions when using this instrument.

Service

In order to guarantee this equipment works safely and efficiently, it must receive regular maintenance. In case of any faults, do not try to repair it yourself. If help is needed, you can always contact your supplier or Labbox via www.labbox.com.

Please provide the customer care representative with the following information:

- Serial number
- Description of problem
- Your contact information

Warranty

This instrument is warranted to be free from defects in materials and workmanship under normal use and service, for a period of 12 months from the date of invoice. The warranty is extended only to the original purchaser. It shall not apply to any product or parts which have been damaged on account of improper installation, improper connections, misuse, accident, or abnormal conditions of operation.

For claim under the warranty please contact your supplier.

Introduction

The V01 mini vortex stirrer adopts advanced international design concepts and manufacturing technology, offering a smooth and attractive appearance, a compact and sturdy structure, high versatility, and simple operability.

This vortex stirrer is used for different samples that require thorough mixing.

Product Overview

Technical Data

Model	V01
Product Type	Vortex stirrer
Control Type	Touch style (fixed speed)
Motor Type	Brushless motor
Orbital Diameter	4.5 mm
Maximum Stirring Capacity	50 ml
Speed Range	3000 rpm
Protection Class	IP42
Output Power	12 W
Voltage	100 – 240 V, 50 / 60 Hz
Permissible Temperature and Humidity	+5 to +40 °C and 80 % RH
Product Size	93 x 93 x 72 mm
Packing Size	190 x 126 x 113 mm
Gross Weight / Net Weight	0.55 kg / 0.36 kg

Operation

- Application: This vortex stirrer is used for different samples that require thorough mixing.
- Usage: It is widely used as an experimental tool in higher education institutions, environmental protection, health, epidemic prevention, the chemical industry, tap water facilities, and other applications.
- Using Steps:
 - Plug in the device and turn the knob to the right to adjust the speed.
 - Place the centrifuge tube or test tube in the top center, then touch and press lightly to start the operation.
-  Instrument Use Environment:
 - Indoor use only
 - Altitude: ≤2000 m
 - Working temperature range: +5 °C to +40 °C
 - Humidity range: ≤80%
 - No vibrations or airflow that could affect performance
 - No conductive dust, explosive gases, or corrosive gases in the surrounding air
 - Maintain a safety distance of 10 cm around the instrument while it is in operation.
-  Safety Precautions:
 - Do not unplug or plug in the power connector or toggle the power button with wet hands.
 - Do not unplug the power plug while the instrument is powered on.
 - Do not perform maintenance or cleaning while the instrument is powered on.
 - Do not install the instrument on an uneven, swaying, or vibrating work surface.

Packing List

Please check the instrument and accessories against the packing list when you first open the packaging. If you find anything wrong with the instrument or accessories, contact the supplier or manufacturer.

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