B10 Vacuum pump

Please read the User Manual carefully before use, and follow all operating and safety instructions!

user manual

english

User Manual

B10 Vacuum pump, LBX B, 10 l/min, 0.075 Mpa

Preface

Users should read this Manual carefully, follow the instructions and procedures, and beware of all the cautions when using this instrument.

Service

In order to guarantee this equipment works safely and efficiently, it must receive regular maintenance. In case of any faults, do not try to repair it yourself. If help is needed, you can always contact your supplier or Labbox via www.labbox.com.

Please provide the customer care representative with the following information:

- Serial number
- Description of problem
- Your contact information

Warranty

This instrument is warranted to be free from defects in materials and workmanship under normal use and service, for a period of 12 months from the date of invoice. The warranty is extended only to the original purchaser. It shall not apply to any product or parts which have been damaged on account of improper installation, improper connections, misuse, accident, or abnormal conditions of operation.


For claim under the warranty please contact your supplier.

Introduction

This B10 vacuum pump adopts advanced international design concepts and manufacturing technology. It features a smooth and attractive shape, a compact and sturdy structure, high versatility, and simple operability, making it a valuable instrument for scientific research experiments.

This vacuum pump is used for vacuum filtrations and other laboratory applications.

Product Overview

Technical Data

| | |
|-------------------------------------|--------------------------------|
| Model | B10 |
| Product Type | Oil-free diaphragm vacuum pump |
| Type of Pressure | Negative pressure |
| Adjustable Pumping Speed | No |
| Pumping Rate | 10 L/min |
| Ultimate Vacuum / Negative Pressure | 0.075 Mpa |
| Noise Level | <40 dB |
| Power | 20 W |
| Material of Valve Diaphragm | Rubber |
| Overload Protection | Yes |
| Overtemperature Protection | Yes |
| Short-Circuit Protection | Yes |
| Product Size | 192 x 100 x 135 mm |
| Weight | 1.58 kg |

Operation

- Application: It is used for vacuum filtrations and other laboratory applications.
- Usage: It is widely used as an experimental tool in higher education institutions, environmental protection, health, epidemic prevention, the chemical industry, tap water facilities, and other applications.
- Using Steps:
 - Plug in the device and turn on the power switch.
-  Instrument Use Environment:
 - Indoor use only
 - Altitude: ≤2000 m
 - Working temperature range: +5 °C to +40 °C
 - Humidity range: ≤80%
 - No vibrations or airflow that could affect performance
 - No conductive dust, explosive gases, or corrosive gases in the surrounding air
 - Maintain a safety distance of 10 cm around the instrument while it is in operation.
-  Safety Precautions:
 - Do not unplug or plug in the power connector or toggle the power button with wet hands.
 - Do not unplug the power plug while the instrument is powered on.
 - Do not perform maintenance or cleaning while the instrument is powered on.
 - Do not install the instrument on an uneven, swaying, or vibrating work surface.

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