

INFORMACIÓN GENERAL

Nombre del producto : Elevador de laboratorio**Descripción :** De aluminio


DATOS TÉCNICOS

referencia	l (mm)	w (mm)	hmáxima (mm)	peso máximo	dimensiones agujero	unidades por ref.
JACA-100-001	100	100	140	30 kg	no tiene	1
JACA-150-001	150	150	280	60 kg	M8x1,25	1
JACA-200-001	200	200	290	60 kg	M8x1,25	1

EMBALAJE Y DATOS LOGÍSTICOS

referencia	vol (l)	kg	TARIC	GTIN
JACA-100-001	1,014	0,477	76169990	08434868027244
JACA-150-001	2,4192	1,27	76169990	08434868027251
JACA-200-001	4,368	1,9528	76169990	08434868027268

FOTO DEL PRODUCTO




GENERAL INFORMATION
Product name : Laboratory jack
Description : Made of aluminum

TECHNICAL DATA

reference	l(mm)	w (mm)	hmax. (mm)	max. weight	hole dimensions	pcs/pack
JACA-100-001	100	100	140	30 kg	does not have	1
JACA-150-001	150	150	280	60 kg	M8x1,25	1
JACA-200-001	200	200	290	60 kg	M8x1,25	1

PACKAGING AND LOGISTICS

reference	vol (l)	kg	TARIC	GTIN
JACA-100-001	1,014	0,477	76169990	08434868027244
JACA-150-001	2,4192	1,27	76169990	08434868027251
JACA-200-001	4,368	1,9528	76169990	08434868027268

PRODUCT PHOTO


INFORMATIONS GÉNÉRALES

Nom produit : Élévateur de laboratoire**Description :** En aluminium

DONNÉES TECHNIQUES

référence	long. (mm)	larg. (mm)	hmaximum (mm)	poids max.	dimensions trou	unités par ref.
JACA-100-001	100	100	140	30 kg	sans	1
JACA-150-001	150	150	280	60 kg	M8x1,25	1
JACA-200-001	200	200	290	60 kg	M8x1,25	1

EMBALLAGE ET LOGISTIQUE

référence	vol (l)	kg	TARIC	GTIN
JACA-100-001	1,014	0,477	76169990	08434868027244
JACA-150-001	2,4192	1,27	76169990	08434868027251
JACA-200-001	4,368	1,9528	76169990	08434868027268

PHOTO PRODUIT



INFORMAZIONE GENERALE

Nome del prodotto : Elevatore regolabile**Descrizione :** In alluminio


DATI TECNICI

referenza	L (mm)	w (mm)	hmassima (mm)	peso massimo	dimensioni foro	unità per ref.
JACA-100-001	100	100	140	30 kg	senza	1
JACA-150-001	150	150	280	60 kg	M8x1,25	1
JACA-200-001	200	200	290	60 kg	M8x1,25	1

IMBALLAGGIO E DATI LOGISTICI

referenza	vol (l)	kg	TARIC	GTIN
JACA-100-001	1,014	0,477	76169990	08434868027244
JACA-150-001	2,4192	1,27	76169990	08434868027251
JACA-200-001	4,368	1,9528	76169990	08434868027268

FOTO DEL PRODOTTO



ALGEMENE INFORMATIE

Produktnaam : Labolift**Beschrijving :** Van aluminium


TECHNISCHE GEGEVENS

referentie	l (mm)	b (mm)	hmaximum (mm)	maximum gewicht	afmetingen opening	stuks per ref.
JACA-100-001	100	100	140	30 kg	niet aanwezig	1
JACA-150-001	150	150	280	60 kg	M8x1,25	1
JACA-200-001	200	200	290	60 kg	M8x1,25	1

VERPAKKING EN LOGISTIEKE GEGEVENS

Referentie	vol (l)	kg	TARIC	GTIN
JACA-100-001	1,014	0,477	76169990	08434868027244
JACA-150-001	2,4192	1,27	76169990	08434868027251
JACA-200-001	4,368	1,9528	76169990	08434868027268

PRODUKTFOTO



ALLGEMEINE INFORMATIONEN
Produktname : Hebebühne
Beschreibung : Aus Aluminium.

TECHNISCHE DATEN

Artikelnummer	L (mm)	B (mm)	Hmaximal (mm)	Gewicht maximal	Lochgröße	Stückzahl pro Artikel
JACA-100-001	100	100	140	30 kg	nicht vorhanden	1
JACA-150-001	150	150	280	60 kg	M8x1,25	1
JACA-200-001	200	200	290	60 kg	M8x1,25	1

VERPACKUNG UND LOGISTIKDATEN

Referenz	vol (l)	kg	TARIC	GTIN
JACA-100-001	1,014	0,477	76169990	08434868027244
JACA-150-001	2,4192	1,27	76169990	08434868027251
JACA-200-001	4,368	1,9528	76169990	08434868027268

PRODUKTFOTO
