

Cell Scrapers

Description

Cell Scrapers, Multi-Position Blade, Non-Pyrogenic, Sterile

Purpose

The process of scraping off and collecting cells more easier and effective

Materials

Handle: ABS (Acrylonitrile-Butadiene-Styrene)


Color: Green

Blade: TPE (Thermoplastic Elastomer)


Color: Natural

Features

- Available with 2 different lengths of 25cm, 39cm
- Two positions of blade: Scraper and Lifter
- These particularly thin, swiveling, flexible blades are easy to use, cells from any damage
- Individually wrapped
- DNase/RNase free and Non-pyrogenic
- Sterilized by irradiation SAL10⁻⁶ (ISO11137)
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 clean room environment
- Manufactured under ISO13485 and ISO9001 quality management system



Dimensions (Unit: mm)



Cell Scrapers

Cat. No.	Total Length	Blade Width	Blade Position	Sterilization	Qty/Pk	Qty/Case
CSCP-025-100	25 cm	2 cm	Scraper	Yes	1	100
CSCP-039-100	39 cm	3 cm	Scraper	Yes	1	100