

Unit: mm

	Dimensions	Tolerance
A	∅4.5	±0.1
B	∅2.45	±0.35
C	268.5	±2.0

SEROLOGICAL PIPETTE

MPIP-J01-500

GRYNIA


DATE: 25-02-2022

Pipette: GPPS (General polystyrene), Compliance with USP Class VI

Filter: PO (Polyolefin)

Colour: Clear

Code colour: Yellow



Unit: mm

	Dimensions	Tolerance
A	∅6.3	±0.1
B	∅2.93	±0.35
C	272.0	±2.0

SEROLOGICAL PIPETTE

MPIP-J02-500

GRYNIA


DATE: 25-02-2022

Pipette: GPPS (General polystyrene), Compliance with USP Class VI

Filter: PO (Polyolefin)

Colour: Clear

Code colour: Green



Unit: mm

	Dimensions	Tolerance
A	∅ 8.2	±0.1
B	∅ 2.5	±0.35
C	341.0	±2.0

SEROLOGICAL PIPETTE

MPIP-J05-500

GRYNIA


DATE: 25-02-2022

Pipette: GPPS (General polystyrene), Compliance with USP Class VI

Filter: PO (Polyolefin)

Colour: Clear

Code colour: Blue



Unit: mm

	Dimensions	Tolerance
A	Ø10.96	±0.10
B	Ø7.95	±0.25
C	Ø2.95	±0.25
D	303.4	±2.0

SEROLOGICAL PIPETTE

MPIP-J10-500

GRYNIA


DATE: 25-02-2022

Pipette: GPPS (General polystyrene), Compliance with USP Class VI

Filter: PO (Polyolefin)

Colour: Clear

Code colour: Orange



Unit: mm

	Dimensions	Tolerance
A	∅15.1	±0.10
B	∅7.95	±0.25
C	∅5.0	±0.25
D	308.5	±2.0

SEROLOGICAL PIPETTE

MPIP-J25-500

GRYNIA


DATE: 25-02-2022

Pipette: GPPS (General polystyrene), Compliance with USP Class VI

Filter: PO (Polyolefin)

Colour: Clear

Code colour: Red



Unit: mm

	Dimensions	Tolerance
A	∅18.4	±0.10
B	∅8.0	±0.25
C	∅5.8	±0.25
D	346.6	±2.0

SEROLOGICAL PIPETTE

MPIP-J50-500

GRYNIA

DATE: 25-02-2022

Pipette: GPPS (General polystyrene), Compliance with USP Class VI

Filter: PO (Polyolefin)

Colour: Clear

Code colour: Purple