

PCR 96-WELL PLATE

**PCRL-096-100**

astik'S


Date: 31-03-2022

Material: Polypropylene USP VI

Colour: Natural, transparent

Capacity: 0,2 ml

Features: Alphanumeric code. Non-skirted. Human DNA free and certified DNase, RNase, and PCR inhibitors free



qPCR 96-WELL PLATE

**PCRL-Q96-100**


astik'S

Date: 31-03-2022

Colour: Natural, transparent

Capacity: 0,2 ml

Features: Alphanumeric code. Semi skirted. Human DNA free and certified DNase, RNase, and PCR inhibitors free



qPCR 96-WELL PLATE

**PCRL-L96-100**

astik'S


Date: 31-03-2022

Material: Polypropylene USP VI

Colour: Natural, transparent

Capacity: 0,1 ml

Features: Alphanumeric code. Sub-semi skirted. Human DNA free and certified DNase, RNase, and PCR inhibitors free



qPCR 96-WELL PLATE

**PCRL-W96-100**

astik'S


Date: 31-03-2022

Material: Polypropylene USP VI

Colour: White, opaque

Capacity: 0,1 ml

Features: Alphanumeric code. Semi skirted. Human DNA free and certified DNase, RNase, and PCR inhibitors free



qPCR 96-WELL PLATE

**PCRL-N96-100**

astik'S


Date: 08-03-2023

Material: Polypropylene USP VI

Colour: Natural, transparent

Capacity: 0,2 ml

Features: Alphanumeric code. Sub-semi skirted. Human DNA free and certified DNase, RNase, and PCR inhibitors free



qPCR 96-WELL PLATE  
**PCRL-P96-100**

astik'S

Date: 08-03-2023

Material: Polypropylene USP VI

Colour: Natural, transparent

Capacity: 0,1 ml

Features: Alphanumeric code. Semi skirted. Human DNA free and certified DNase, RNase, and PCR inhibitors free