

Slide mailers (PP material)

1. Technical information		
Product name	Slide mailers (PP material)	
Product description	Slide mailers, PP material, available for 1, 2, 3, 5 standard 3 x 1" slides, available in edge opening or end opening. Reusable.	
2. Intended Use		
The slide mailer made of PP material is ideal for shipping, storing and handling of microscope slides.		
3. Features		
<ul style="list-style-type: none"> • Made from PP material, break-resistance. • Handy, convenient and economical • Hold microscope slides in place so they will not touch when the cover is closed. • Tight seal keeps away from dust and moisture. • Supplied with Neutral and Blue color. 		
4. Certification		
It is certified that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our EN ISO13485:2016&EN ISO9001:2015 Quality Management System.		
5. Specification		
P/N	Specification	Packaging(box/carton)
MSVP-005-100	5-pcs, Neutral color.	100pcs/dispenser box, 1,000pcs/carton
6. Sterilization		
<input type="checkbox"/> STERILE EO <input type="checkbox"/> STERILE R <input type="checkbox"/> STERILE A <input checked="" type="checkbox"/> Non-sterile		
7. Precautions		
<ul style="list-style-type: none"> • Clean up after use if necessary reused. 		
8. Storage and Shelf Life		
Store in a clean, dry environment at room temperature.		