

COSS-10U-500**METRIA**

Temperature °C	Conductivity uS/cm	Temperature °C	Conductivity uS/cm
5	6.13	22	9.38
10	7.10	23	9.59
15	7.95	24	9.79
16	8.15	25	10.00
17	8.37	30	11.03
18	8.56	35	12.14
19	8.78	40	13.29
20	8.97	45	14.44
21	9.18	50	15.55

COSS-084-500

Temperature °C	Conductivity uS/cm	Temperature °C	Conductivity uS/cm
5	50.40	22	78.96
10	58.80	23	80.64
15	67.20	24	82.32
16	68.88	25	84.00
17	70.56	30	92.40
18	72.24	35	100.8
19	73.92	40	109.2
20	75.60	45	117.6
21	77.28	50	126.0

COSS-147-500

Temperature °C	Conductivity uS/cm	Temperature °C	Conductivity uS/cm
5	88.20	22	138.18
10	102.90	23	141.12
15	117.60	24	144.06
16	120.54	25	147.00
17	123.48	30	161.70
18	126.42	35	176.40
19	129.36	40	191.10
20	132.30	45	205.80
21	135.24	50	220.50

COSS-1K4-500**METRIA**

Temperature °C	Conductivity uS/cm	Temperature °C	Conductivity uS/cm
5	894	22	1331
10	1007	23	1358
15	1139	24	1386
16	1167	25	1413
17	1194	30	1550
18	1221	35	1694
19	1249	40	1833
20	1276	45	1989
21	1304	50	2139

COSS-012-500

Temperature °C	Conductivity uS/cm	Temperature °C	Conductivity uS/cm
5	7728.0	22	12107.2
10	9016.0	23	12364.8
15	10304.0	24	12622.4
16	10561.6	25	12880.0
17	10819.2	30	14168.0
18	11076.8	35	15456.0
19	11334.4	40	16744.0
20	11592.0	45	18032.0
21	11849.6	50	19320.0

COSS-111-500

Temperature °C	Conductivity uS/cm	Temperature °C	Conductivity uS/cm
5	67080	22	105092
10	78260	23	107328
15	89440	24	109564
16	91676	25	111800
17	93912	30	122980
18	96148	35	134160
19	98384	40	145340
20	100620	45	156520
21	102856	50	167700